

CLINICON®

Laboratory Equipment & Supplies

Pluggo Abbot Alinity - Technical description

Pluggo is designed for safe removal of original caps from vacuum collection tubes. Pluggo is capable of handling 5ml, 7ml, 7ml long and 10ml vacuum tubes by use of the three different coloured carousels.

PLUGGO



The caps are guided to an outlet for collection of biohazard disposal with a detachable chute.

Technical data

Part number:	9600 7450
Sample tubes:	Vacuum collection tubes
Dimensions including carousel:	430x270x210 mm (Depth x Width x Height)
Weight:	9.5 Kg
Voltage:	230VAC 50Hz / 110VAC 60Hz 60VA