


Pluggo Standard Automated Decapper - Technical description

Pluggo is designed for safe removal of original caps from vacuum collection tubes. Pluggo is capable of handling 5ml, 7ml, 7ml long and 10ml vacuum tubes by use of the three different coloured carousels.

Each carousel is capable of handling up to 50 sample tubes.

PLUGGO



Carousel colour	Tube size	Volume	Number of tubes
Blue	13 x 75 mm	5 ml	50
Red	13 x 100 mm	7 ml long	50
Yellow	16 x 100 mm	10 ml	50

The caps are guided to an outlet for collection of biohazard disposal with a detachable chute.

Technical data

Part number:	9600 7950
Capacity:	Up to 50 samples in each carousel
Sample tubes:	Vacuum collection tubes
Speed:	55 seconds per carousel
Dimensions including carousel:	580x350x210 mm (Depth x Width x Height)
Carousel dimension (ex handle)	90x350 mm (Height x Diameter)
Weight:	9.5 Kg
Weight of carousel	1.6 Kg
Voltage:	230VAC 50Hz / 110VAC 60Hz
Power:	60VA