

Technical Data Sheet

Supplier:	Clinicon AB
Product Name:	Neutral over-cap for 13mm tube
Product Number:	2174-3211
Product Description:	Push-on over cap for laboratory use
Material of the cap:	NOVEX LM2010AA LDPE polyethylene
Storage Condition (°C):	N/A
Expiry:	N/A
Sanitary Quality:	Clean
Scope of uses:	Industries, Clinical
Unit of Sale:	1000 (BULK PACKED IN 20,000)
Net weight:	340g/1000
Gross weight:	6.83kg(BAG20,000)

Clinicon products are sold through Clinicon AB. Clinicon warrants that its products conform to the information contained in this and other Clinicon publications. Purchaser must determine the suitability of the product for their use.

Regulation Statement

We hereby confirm that we do not intentionally use the following substances in the manufacturing of the product:

-LATEX

-BISPHENOL A (CAS No 80-05-7)

-PHTHALATES as per Regulation 1907/2006/EC (SVHC)

Raw material suppliers also guarantee the following requirements through certificates of compliance:

-EU Food Contact 1935/2010/EC and 10/2011/EC

We hereby confirm that we do not intentionally use the following substances in the manufacturing of the product: prioniform encephalopathy agents via human and veterinary medicinal products* by virtue of all bovine derived material having

been processed per specific conditions of section 6.4 of EMA/410/01.

-Compliance to REACH SVHC candidate list (Regulation 1907/2006/EC). They certify that the here above mentioned products do

not contain one than 0.1% (w/w) of any substance of very high concern listed in the ECHA candidate list in accordance with

Article 33 of the Reach regulation as updated on January 12, 2017.

-Heavy Metals Free as per directive 94/62/EC. The chemical composition of this product meets the relevant requirements of

Directive 92/62/EC (as amended) on packaging and waste. The sum of the concentration Cd, Cr (V1), Hg and Pb in this product

does not exceed 100 ppm.

Raw material suppliers also confirm through certificates of compliance that they do not intentionally use the following substances:

-LATEX

-BISPHENOL A (CAS No 80-05-7)

-PHTHALATES as per Regulation 1907/2006/EC (SVHC)

Clinicon products are sold through Clinicon AB. Clinicon warrants that its products conform to the information contained in this and other Clinicon publications. Purchaser must determine the suitability of the product for their use.