


Pluggo Standard Automated Decapper - Technical description

Pluggo is designed for safe removal of original caps from vacuum collection tubes. Pluggo is capable of handling 5ml, 7ml, 7ml long and 10ml vacuum tubes by use of the different coloured carousels.



PLUGGO



Carousel colour	Tube size	Volume
Blue	13 x 75 mm	5 ml
Red	13 x 100 mm	7 ml long
Grey	16 x 75 mm	7 ml short
Yellow	16 x 100 mm	10 ml

The caps are guided to an outlet for collection of biohazard disposal with a detachable chute.

Technical data

Capacity:	Up to 24 samples in each carousel
Speed:	30 seconds per carousel
Size:	430x270x210 mm (Depth x Width x Height)
Weight:	9,0 kg
Voltage:	230VAC 50Hz / 110VAC 60Hz
Power:	60VA